

Recommended use of heat + electrodes

The appropriate selection of ECG electrodes depends on the test type, the patients' age and their physical features. The variety of materials and gels used for the production of SORIMEX heart+ electrodes allows the use of appropriate electrodes for the individual needs of medical personnel and patients.

type of medical examination	EK-1/23 WpHG/2j-L50	EK-1/23 WpHG/2j-L50/RTG	22 PSG	25 PSG	30 PSG	35 PSG ST	40 PSG	42 PSG	44 PSG	45 PSG	45 PSG RTG	50 P	50 PSG	50 NWSG	55 P	55 PSG	56 P	60 P	60 PSG	60 PSG XC	60 PT	61 P	61 PSG
long-term	x	x									x			x				x	x	x	x	x	x
short-term	x	x	x	x	x			x	x		x												
resting			x	x	x		x	x	x	x		x	x		x	x	x						
monitoring	x		x	x	x		x	x	x	x		x	x		x	x	x						
diagnostic	x		x	x	x		x	x	x	x		x	x		x	x	x						
RTG examination		x									x												
Holter examination														x				x	x	x	x	x	x
Stress-test						x				x					x								
examination at ICU														x									
continous cardiac supervision														x									

